

Technically
Speaking

CEDARLANE[®] 
www.cedarlanelabs.com

Conveniently Delivering You Today's Innovations
for the Science of Tomorrow™

FLAER: Alexa fluor[®] 488 proaerolysin

FL1S – 25 µg FLAER in 0.5 ml PBS-albumin

FL2S – 50 µg FLAER in 1.0 ml PBS-albumin

Description:

The concentration of FLAER in the supplied stock solution is 10^{-6} M. The solution also contains 5 mg/mL albumin, added in order to maximize stability. **The stock solution may be kept in the refrigerator protected from light for at least 1 year.**

FLAER is produced by labeling an inactive variant of the protein proaerolysin with Alexafluor[®] 488 succinimidyl ester, provided under agreement with Life Technologies Corporation. The reagent binds selectively to cell surface GPI-anchored proteins. It is intended for research purposes only and is not for resale.

5×10^{-8} M has proven to be a useful final concentration of FLAER for flow cytometry (this represents a 1 to 20 dilution of the stock directly into the test sample), however you may wish to try other concentrations for your own application.

There are a number of published methods describing the detection of PNH using FLAER. One example is Sutherland et al. in Cytometry Part B (Clinical Cytometry) 72B: 167-177 (2007). We can provide a detailed protocol upon request.

Manufactured by:

PINEWOOD SCIENTIFIC SERVICES INC
210 Denison Rd., Victoria, BC, V8S 4K3, Canada
Tel: 250-598-6822, Fax: 250-598-6877
protoxbiotech.com

Updated: January 17, 2011

Visit our website for your local distributor.

CEDARLANE[®] 

www.cedarlanelabs.com

An ISO 9001:2000 and ISO 13485:2003
registered company.

In CANADA: **Toll Free: 1-800-268-5058**

4410 Paletta Court, Burlington, ON L7L 5R2 ph: (289) 288-0001, fax: (289) 288-0020
e-mail: general@cedarlanelabs.com

In the USA: **Toll Free: 1-800-721-1644**

1210 Turrentine Street, Burlington, NC 27215 ph: (336) 513-5135, fax: (336) 513-5138
e-mail: service@cedarlanelabs.com