

TETANUS TOXOID

FOR IN VITRO TESTS



Highly purified Tetanus Toxin from *Clostridium tetani*, which following the purification process has been detoxified by formaldehyde.

Product Description

- The product is supplied with a toxoid concentration of 500 – 1.200 Lf/ml, corresponding to a specific protein concentration of 1.25 – 3.00 mg/ml.
- The product can be regarded as close to a 100% pure.
- The buffer is saline.
- The concentration of formaldehyde is ≤ 75 mg/l.
- The product is supplied in 10 ml capped vials.
- The recommended storage temperature is 2-8°C.
- The product is shipped at ambient temperature.

Applications

The product is suitable for various purposes, e.g.:

- Lymphocyte stimulation tests (proliferation tests).
- Induction and measurements of cytokine synthesis.
- Studies of antigen induced migration of monocytes/macrophages and other phagocytizing cells.
- ELISA

Art No.	Item
2674 (SS)	Tetanus Toxoid for in vitro, 10 ml
101778 (SS)	Tetanus Toxoid for in vitro, 1 ml

These products used to be sold by Statens Serum Institut



Part of TCP Analytical Group

CANADA

4410 Paletta Court
Burlington, ON, Canada, L7L 5R2
E-mail: general@cedarlanelabs.com
Toll Free: 1-800-268-5058
Ph: (289) 288-0001
Fax: (289) 288-0020

USA

1210 Turrentine St, Burlington, North Carolina - 27215, United States
E-mail: service@cedarlanelabs.com
Toll Free: 1-800-721-1644
Ph: (336) 513-5135
Fax: (336) 513-5138