

For research use only

AcidiFluor ORANGE™-Dextran 10k

Code No.	Material	Amount	Storage	Stability
GC306	AcidiFluor ORANGE™- Dextran 10k	1 mg	Store at -20 °C, desiccate and protect from light	1 year (unopened)

AcidiFluor ORANGE™-Dextran 10k is prepared by conjugating AcidiFluor ORANGE™ to aminodextran. This probe has weak fluorescence at neutral pH, but shows strong fluorescence in acidic conditions such as in acidic organelles, such as, lysosomes and late endosomes. It can be used to detect endocytosis.