

For research use only

# AcidiFluor ORANGE™-Zymosan A

---

Code No.	Material	Amount	Storage	Stability
GC305	AcidiFluor ORANGE™- Zimosan A	1 mg × 1 vials	Store at -20 °C, desiccate and protected from light	1 year (unopened)

**AcidiFluor ORANGE™-Zymosan A** is prepared by conjugation of AcidiFluor ORANGE™ to Zymosan A. This particle shows weak fluorescence at neutral pH, but shows strong fluorescence in acidic organelles such as lysosomes and late endosomes. It can be used to detect phagocytosis.