

## ***Diaminofluorescein-4M acetoxymethyl ester (DAR-4M AM)***

Catalog no.	Product name	Amount	Storage upon receipt	Stability
SK1006-01	Diaminorhodamine-4M acetoxymethyl ester (DAR-4M AM)	1 mg (in DMSO 0.32 ml)	Store frozen, desiccate and protect from light.	1 year (unopened)

**DAR-4M AM** is an acetoxymethyl ester derivative of diaminorhodamine-4M. Cell permeable DAR-4M AM is hydrolyzed by intracellular esterase to generate cell impermeable DAR-4M, thus, it can detect intracellular NO in pH range of 4 – 12.

