

For research use only

## STELLA Fluor™720-NHS

Catalog no.	Product name	Amount	Storage upon receipt	Stability
ST1011-11 ST1011-15	STELLA Fluor™ 720-NHS	1 mg 5 mg	Store at <-20°C, desiccate and protect from light.	1 year (when unopened and stored as described.)

**STELLA Fluor™720-NHS** is an amine-reactive fluorescent dye. It can be used to label proteins or other amine-containing chemicals. By mixing with the proteins, it forms stable covalent bond with lysine residue. A condensation reagent is NOT required for the reaction. 1 mg of this product will label  $2.1 - 6.3 \times 10^{-7}$  mol of protein.