

Technically
Speaking

CEDARLANE[®]
www.cedarlanelabs.com



Conveniently Delivering You Today's Innovations
for the Science of Tomorrow™

**Anti-Human IgG Fc
Monoclonal Antibody**

Catalogue#	Format	Size	Concentration	Isotype Control
CL6017AP	Purified	50 µg	1.0 mg/ml	CLCMG2A00
CL6017AP-2	Purified	250 µg	1.0 mg/ml	CLCMG2A00
CL6017AP-3	Purified	500 µg	1.0 mg/ml	CLCMG2A00
CL6017B	Biotin	100 µg	0.1 mg/ml	CLCMG2A15
CL6017F	FITC	100 µg	0.1 mg/ml	CLCMG2A01
CL6017PE	PE	50 µg	0.1 mg/ml	CLCMG2A04
CL6017HP	HRPO	100 µg	1.0 mg/ml	CLCMG2A07

Isotype: Mouse IgG_{2a}

DESCRIPTION:

The Fc region of IgG is a homodimer composed of the constant region of the two heavy chains that form the IgG molecule. The Fc region mediates a number of processes including opsonization, antibody dependent cellular cytotoxicity (ADCC), and complement activation through binding to Fc receptors such as CD16, CD32, CD64, and the complement factor C1.

PRESENTATION:

Purified: Purified IgG buffered in PBS and 0.02% NaN₃. (Purified from ascitic fluid via Protein G Chromatography). For maximum recovery of contents, spin down tube before use.

Biotin, FITC, and PE: Biotin/FITC/PE conjugated IgG buffered in PBS, 0.02% NaN₃ and EIA grade BSA as a stabilizing protein to bring total protein concentration to 4-5 mg/ml.

HRPO: HRPO conjugated IgG buffered in PBS with 40% glycerol (v/v) and EIA grade BSA as a stabilizing protein to bring total protein concentration to 4-5 mg/ml. No NaN₃ with other preservatives.

STORAGE/STABILITY:

For all formats, store at + 4°C. DO NOT FREEZE PE conjugates. For long term storage (**Purified/Biotin/FITC/HRPO**), aliquot and freeze unused portion at -20°C in volumes appropriate for single usage. Avoid freeze/thaw cycles. For Maximal recovery of contents, please quick-spin vial before opening.

Continued Overleaf.....

Visit our website for your local distributor.

CEDARLANE[®]



www.cedarlanelabs.com

An ISO 9001:2000 and ISO 13485:2003
registered company.

In CANADA: Toll Free: 1-800-268-5058

4410 Paletta Court, Burlington, ON L7L 5R2 ph: (289) 288-0001, fax: (289) 288-0020
e-mail: general@cedarlanelabs.com

In the USA: Toll Free: 1-800-721-1644

1210 Turrentine Street, Burlington, NC 27215 ph: (336) 513-5135, fax: (336) 513-5138
e-mail: service@cedarlanelabs.com

APPLICATION:

This antibody is suitable for use in ELISA (both capture and detection antibody), Radioimmunoassay, Western Blot and Flow Cytometry.

*** For optimal results in various applications, it is recommended that each investigator determine dilutions appropriate for individual use.**

SPECIFICATIONS:

Clone: HP6017

Immunogen: Purified human immunoglobulin

Specificity: Human IgG Fc

ELISA Specificity: human IgG subclass 1 100%
human IgG subclass 2 100%
human IgG subclass 3 100%
human IgG subclass 4 99%
human IgG Pan (pooled IgG1-4) 100%
human IgG Fc 100%
human IgG Fab 0.12%

REFERENCES:

1. Hamilton, R.G. et al; Isoelectric focusing-affinity immunoblot analysis of mouse monoclonal antibodies to the four human IgG subclasses. Electrophoresis, 8:127-34, 1987.
2. Papadea, C. et al; Monoclonal antibody-based solid-phase immunoenzymometric assays for quantifying human immunoglobulin G and its subclasses in serum. Clin. Chem., 31: 1940-5, 1985.
3. Jeffries, R. et al; Evaluation of monoclonal antibodies having specificity for human IgG sub-classes: results of an IUIS/WHO collaborative study; Immunol Lett. 10(3-4):223-52, 1985.
4. Reimer, C.B. et al; Evaluation of thirty-one mouse monoclonal antibodies to human IgG epitopes. Hybridoma, 3: 263-75, 1984.

FOR RESEARCH USE ONLY