


Urinary Creatinine kits N-Cal™ NIST-Calibrated Catalog No. K002-H1/H5

KIT FEATURES


- US NIST-Calibrated Standard
- Safe, No-Mix Detection Reagent
- Validated in a Variety of Samples
- 30 Minute Assay, 4°C Stable Reagents
- 88 or 472 samples in duplicate



Urea Nitrogen (BUN) kit N-Cal™ NIST-Calibrated Catalog No. K024-H1

KIT FEATURES


- US NIST-Calibrated Standard
- Sensitivity = 30 µg/dL
- 30 Minute Assay, 4°C Stable Reagents
- 88 samples in duplicate



Dual Range Hemoglobin kit Catalog No. K013-H1

KIT FEATURES

- Blood, RBCs, Serum, Plasma
- Sensitivity = 20 µg/dL
- 30 Minute Assay
- 88 samples in duplicate




Phone: 734.677.1774
E Mail: Info@ArborAssays.com
www.ArborAssays.com



Cortisol EIA kits N-Cal™ NIST-Calibrated Catalog No. K003-H1/H5

KIT FEATURES


- Urine, Fecal Extracts, Saliva, plus:
- 1.5 µL Serum/Plasma Sample for Duplicate Measurements
- Tested in Human, Ruminants, Fish, Cats, Rodents, Birds
- 90 Minute Assay, Stable Reagents
- 39 or 231 samples in duplicate



Corticosterone EIA kits Catalog No. K014-H1/H5

KIT FEATURES

- Fecal Extracts, Saliva, TCM plus:
- 1.5 µL Serum/Plasma Sample for Duplicate Measurements
- Tested in Human, Ruminants, Fish, Cats, Rodents, Birds
- 90 Minute Assay, Stable Reagents
- 39 or 231 samples in duplicate




Phone: 734.677.1774
E Mail: Info@ArborAssays.com
www.ArborAssays.com



Cortisone Chemiluminescent (CLIA) kits Catalog No. K017-C1/C5

KIT FEATURES


- Measure Activity of 11 β -HSD isozymes
- 1.5 μ L Serum Sample for Duplicate Measurements
- Tested in Human, Ruminants, Fish, Cats, Rodents, Birds
- 120 Minute Assay, 4 $^{\circ}$ C Stable Reagents
- 37 or 229 samples in duplicate
- World's Only Cortisone Immunoassay



13,14-dihydro-15-keto-PGF_{2 α} (PGFM) EIA kits Catalog No. K022-H1/H5

KIT FEATURES


- Pregnancy Marker in Felids
- Non-Invasive Urine and Fecal Extracts
- Plus Serum, Plasma and TCM
- 90 Minute Assay, 4 $^{\circ}$ C Stable Reagents
- 39 or 231 samples in duplicate



Progesterone EIA kits Catalog No. K025-H1/H5

KIT FEATURES


- Measures most 3- and 11-hydroxy Progesterone metabolites
- Urine, Fecal Extracts, TCM
- Tested in Human, Ruminants, Fish, Cats, Rodents, Birds
- 150 Minute Assay, Stable Reagents
- 39 or 231 samples in duplicate



17 β -Estradiol EIA kits Catalog No. K030-H1/H5

KIT FEATURES

- Uses the Most Published Antibody for Urinary and Fecal E2 Detection
- Tested in Human, Ruminants, Fish, Cats, Rodents, Birds
- 150 Minute Assay, Stable Reagents
- 41 or 233 samples in duplicate



Phone: 734.677.1774
 E Mail: Info@ArborAssays.com
www.ArborAssays.com



Phone: 734.677.1774
 E Mail: Info@ArborAssays.com
www.ArborAssays.com

